



Quanta System

◊ SUPREMA ◊

AESTHETICS & DERMATOLOGY

Quanta System

◇ SUPREMA ◇

UNLEASH THE SUPREMA POWER

Every pulse is like the first light of the day, designed to **renew the skin's innate radiance and vitality.**

Suprema is for those who are seeking treatments that reveal their **authentic beauty.**

BOOST YOUR CHANCES

Suprema is designed to **meet multiple skin needs**, allowing clients to experience a profound **sense of rejuvenation.** This **commitment to well-being** has enabled a wide range of aesthetic procedures.

ADVANCED TECHNOLOGY

This laser device has two main features at its core: **three selectable wavelengths** and the **Quanta Mixed Technology®**, which allows **two different wavelengths** to be combined with **large spot sizes.**

RAISE YOUR COMFORT

Suprema's ability to address a **variety of different procedures**, along with the presence of a comfort enhancing double **Cooling System**, allows consumers to find a solution for their own skin needs.

OVERSTEP YOUR LIMITS

A **versatile laser device** like Quanta System's Suprema is the perfect addition to your clinic if you want to **elevate the overall quality of treatments** and **increase your return on investment.**

Overview

Up to 3 Wavelengths

Nd:YAG 532 nm

Nd:YAG 1064 nm

Alexandrite 755 nm

Quanta Mixed Technology®

Proprietary technology that allows to mix and emit wavelengths Nd:YAG (1064 nm) and Alexandrite (755 nm) simultaneously or sequentially

Contact Skin Cooler Handpiece

Specific handpiece that enables the use of contact cooling, enhancing the comfort sensation during treatment

Air Skin Cryo Handpiece

Innovative handpiece designed to enhance comfort working in combination with an external air cooling unit

Optional Twain Handpieces

Twain Er:YAG 2940 nm and Twain IPL can be connected to further extend the treatment possibilities



A POWERFUL SOLUTION FOR MULTIPLE APPLICATIONS

A woman with long, wavy brown hair is shown from the chest up, wearing a purple bikini top. She is looking upwards and to the right, with her eyes closed, against a background of a calm ocean and a hazy sky. The text "Boost Your Chances" is overlaid on the left side of the image.


Boost Your
Chances

Versatility is the — real key

At Quanta System, we are dedicated to raising awareness about the importance of **embracing diversity and celebrating the unique beauty of every individual**.

Suprema was designed for **high versatility**, addressing multiple skin needs such as **dermatological vascular lesions, benign pigmented lesions, hair removal** and **skin rejuvenation** through specific procedures.

This **advanced laser device** unlocks new possibilities in aesthetic treatments, thanks to its **remarkable key features** and the **innovative technology** that distinguish it in the contemporary laser market.



VARIETY OF APPLICATIONS

DERMATOLOGICAL
VASCULAR
LESIONS

HAIR
REMOVAL

BENIGN
PIGMENTED
LESIONS

SKIN
REJUVENATION

TAKING CARE OF EVERY SKIN NEED

Find the right — combination

Suprema is a **smart choice for laser aesthetic treatments** thanks to its advanced features. Specifically, this powerful laser device leverages the full power of **Quanta Mixed Technology®**.

This technology allows **the combination of two different laser wavelengths** in a single pulse, providing **versatility and effectiveness** across various treatments and conditions.

With **Quanta Mixed Technology®**, Suprema specifically enables the **simultaneous or sequential mixing and emission of two wavelengths: Alexandrite (755 nm) and Nd:YAG (1064 nm)**.


Through this **meticulous customization**, procedures can be **tailored to the specific needs of the consumer**, enhancing treatment safety and efficiency, and thereby delivering **high-level outcomes**.



SUPREMA
Quanta Mixed Technology®

ALEXANDRITE (755 nm)

ND:YAG (1064 nm)



THE RIGHT FORMULA FOR THE BEST PROCEDURES



Advanced Technology

Raise your — comfort

With Suprema, we deeply focus on **wellness preservation** in every procedure, providing efficient solutions that can truly **transform consumers' perception of their treatments**.

In this journey, **technological innovation** plays an important role, enabling Suprema to **address a wide variety of skin concerns**, significantly understanding people's needs, and **always being on the customer's side**.

This caring approach provides consumers with **a sense of regeneration**, allowing them to **experience comfort in an authentic way**, leading to **a profound sense of well-being**.

This is your moment to feel good.

Unleash your Suprema power.

Dermatological Vascular Lesions

Hair Removal

Benign Pigmented Lesions

Skin Rejuvenation



Treatments



Dermatological Vascular Lesions

By leveraging advanced laser technology, **damaged blood vessels** get selectively targeted and closed, reducing vascular lesions while **preserving surrounding tissues**.

The **laser's interaction with hemoglobin** in the capillary vessels effectively leads to their coagulation and subsequent disappearance.

Dermatological vascular lesions that can be treated include **telangiectasias, varicose veins, hemangiomas, rosacea** and **vascular malformations**, among others.




Hair Removal

Facial and body unwanted hair can be certainly annoying in many situations, both for aesthetic and dermatological reasons.

Laser technology addresses **hair removal necessities** by targeting hair follicles and acting on **melanin concentration**, while **minimizing damage** to surrounding skin.

During the treatment, **hair follicles** absorb the pulses of light emitted by the device. This procedure removes the hair, **smoothing skin** and **inhibiting future hair growth**.



Benign Pigmented Lesions

Benign pigmented lesions can be caused by a **variety of factors** and certainly need to be treated, even if they are not dangerous to consumer's health.

These procedures efficiently address **targeted skin areas**, meeting **customers' specific needs**.

Well-being is preserved and efficacy is assured on a variety of indications among which **melasma, lentigos, freckles, café au-lait spots and Ota & Ito nevus**.



Skin Rejuvenation

While aging is a natural process for everyone, laser skin rejuvenation offers a real solution, targeting various skin concerns and **specific areas of the face and neck**.

This procedure enables timely action by **inducing collagen production** and **reducing a wide range of aging signs**.

Skin rejuvenation is effective for treatment on **all skin types**, providing consumers with a **smoothing effect** as a result of **multi-level skin stimulation**.

Handpieces



Contact Skin Cooler

Contact cooling for comfort enhancement



Interchangeable Zoom Optics

Two options available for addressing diverse procedures

Available Dimensions

2,5-6 mm

8-20 mm



Twain 2940

Er:YAG laser @ 2940 nm. The ablative power in your hands

Air Skin Cryo

The full power of air cooling in an innovative handpiece



Interchangeable Fixed Optics

Five different spot dimensions for highly accurate results

Available Spot Sizes

14 mm
16 mm
18 mm
20 mm
22 mm



Twain IPL

The mix of power and versatility of Intense Pulsed Light





Overstep
your Limits

Make your — best choice

Suprema has been designed by Quanta System to **enhance aesthetic and dermatological procedures**, enabling operators to perform **a wide range of treatments**.

With its notable potential, Suprema can **make a significant difference in the aesthetic market**, helping clinics and operators **build strong reliability** while meticulously **caring for consumers**.

Quanta Mixed Technology[®] is certainly the key feature you need if your ambition is to overstep the limits and **raise the bar in aesthetic procedures, upgrading your own practice**.



*PULSES
COMBINATION*

*HIGH LEVEL
PROCEDURES*

*OUTCOMES
IMPROVEMENT*

EXPLORE NEW TREATMENT POSSIBILITIES

Before/After



Dermatological Vascular Lesions Rosacea

Courtesy of Nicola Zerbinati MD



Dermatological Vascular Lesions Telangiectasia

Courtesy of Nicola Zerbinati MD



Dermatological Vascular Lesions Hemangioma

Courtesy of Nicola Zerbinati MD



Hair Removal Beard Design

Courtesy of Quanta System



Benign Pigmented Lesions Solar Lentigo

Courtesy of Nicola Zerbinati MD



Advantages

COMBINED TECHNOLOGY

With Quanta Mixed Technology[®], operators can reach top-level results through the powerful combination of two wavelengths

COMPETITIVE COVERAGE

Suprema enables access to a wide range of treatment possibilities, increasing the clinic's competitive power

CLIENT-ORIENTED SOLUTIONS

Thanks to the double Cooling System, Suprema cares for customers' well-being while preserving their sensation of comfort

IMPROVED REVENUES

Suprema has the potential to be the top-choice device you are looking for to significantly upgrade your practice

MORE INFO HERE



Quanta System

get in touch



[quantasystem.com](https://www.quantasystem.com)



This brochure is not intended for the U.S. market.
Certain Intended Uses/Configurations/Models/Accessories are not cleared for U.S.

DANGER - VISIBLE AND INVISIBLE LASER RADIATION

Avoid eye or skin exposure to direct or scattered radiation | Laser product: Class 4



Quanta System S.p.A. – Via Acquedotto, 109, 21017 Samarate (VA) – Italy
© All rights reserved – quanta@quantasystem.com | www.quantasystem.com